



Zetta Medical Technologies
 1313 Ensell Rd.
 Lake Zurich, IL 60047
 Tel. 847-550-9990
 Fax. 847-550-9994
<http://www.zettamed.com>

Philips MX 8000 Dual CT

- Two (2) systems available
- **Manufactured 2001**
- Water Chiller included
- Computer UPS Included
- Extremely Clean - like new condition
- Available Immediately
- Delivery available worldwide
- Installation and calibration
- Local service delivery
- 12 months parts, tube and labor warranty included
- Systems fully tested and calibrated before shipping.



Installation & Warranty:

Zetta will install and calibrate the two systems. 12 months warranty for labor, parts and CT tube is included.

Service Coverage:

- Labor & travel from 8:00AM – 5:00PM Monday through Friday (Excluding Holidays)
- Priority Onsite Service Response
- Four (4) scheduled Planned Maintenance visits per year per system
- Replacement parts are included
- Priority Parts Delivery
- Detailed Service Management Reports

Model	Mx8000D	Mx8000D
System ID	ZM5050	ZM5121
Generator	48KW	48KW
Specifications	Details	
Scan Time	0.75, 1.0, 1.5, 2.0 sec	
Slice Width	2 X 1, 2 X 2.5, 2 X 5, 2 X 8, 2 X 10	
Gantry Tilt	+/- 30 Deg	
Resolution	Standard, High, Ultrahigh	
FOV	500 mm	
Gantry Aperture	700 mm	
Patient Support	200kg/ 440 Lbs	
DICOM	Supports SCU, SCP, DICOM Print, Modality Worklist, Query/Retrieve,	
Computer	Silicon Graphics O2	
Software Options	Host Computer	Host Computer
Viewer	Yes	Yes
MPR	Yes	Yes
3D	Yes	Yes
Angio-MIP	Yes	Yes
DentaCT	No	No
MasterMatch	No	No
Time Lapse	Yes	Yes
Combine Images	Yes	Yes
Voyager	No	No
4-D Angio	Yes	No
Heartbeat CS	No	Yes
Functional CT	No	No
QCTA	Yes	Yes
MasterLook	Yes	Yes
Stereotaxis	Yes	Yes
MRI Normalisation	Yes	Yes