



Zetta Medical Technologies, LLC
1313 Ensell Rd., Lake Zurich, IL 60047
Tel. 847.550.9990 Fax 847.550.9994
www.ZettaMed.com

2003 Siemens Sensation 16 CT Scanner

16 Slice

Equipment Description:

Recently de-installed, 2003 SIEMENS Sensation 16 CT

- DURA Akron Q High Performance 5.3 MHU X-Ray Tube
- Ultra Fast Ceramic (UFC) Detector
 - Acquiring up to 16 slices per rotation
- 60 kW High Frequency X-Ray Generator
- WorkStream Computer System
- Wizard Console
- (2) 18" High Resolution Color LCD Monitors
- with a 1280 X 1024 Display Matrix
- Operating Syngo Software (CT 2006G)
- Patient Registration, Examination Card,
- Viewing Card, Filming Card
- CARE Solutions:
 - CARE Filter, CARE Topogram, CARE
 - Bolus, CARE Dose, Pediatric Scan Protocols
- Image Acquisition Modes:
 - Topogram, Sequential Scan, Spiral Scan, Dynamic Multi-Scan, Serio-Scan
- Image Display & Evaluation Mode:
 - Image Reconstruction, Image Display, Cine Display, Multiplanar Reformatting, 3D, 3d
- Display, CT Angio
 - DICOM (Query/Retrieve, Storage, Print, Worklist)
- HeartView
- Water Chiller



This system is available immediately and it can be delivered and installed at your facility with a comprehensive 12 months warranty and extended service coverage

Accessories, Injectors, UPS and Printers are also available upon request. Remember to ask about our PACS systems.