



Zetta Medical Technologies, LLC
1313 Ensell Rd., Lake Zurich, IL 60047
Tel. 847.550.9990 Fax 847.550.9994
www.ZettaMed.com

2008 Toshiba Titan 1.5T MRI System

16 Channel System

Software Version V9.50*R241

Equipment Description:

Includes 71 cm bore, 55x55x50 cm FOV, Pianissimo Noise Reduction Technology, Toshiba's proprietary contrast-free MRA techniques (FBI, CIA, Time-SLIP)

Software Licenses Included:

AMB	Noise Reduction
AutoVoice	PPS_SCU
BIP	PRINTSCU
COMMTSCU	QUERYSCP
DBDT02_1	QUERYSCU
DICOM	SAR02_01
DICOMDIR	SILENT
EPI	SPEEDER
JET	STRG_SCU
MRA	SUPER_FASE
MWM_SCU	

Coils Included:

Atlas Speeder Body Coil	Atlas Speeder Body Coil (ext)
Atlas Speeder Spine Coil	Rectangular Flex Coil
Head Support Assembly	Head Coil Bridge
100 Flex Coil	4-Channel Array Adaptor
Head Coil Bridge	Speeder Breast Coil
QD Head Coil	Atlas Speeder Head/Cervical Coil
200 Flex Coil	Wrist Speeder Coil
Shoulder Array Coil	QD Knee/Foot Coil



This system is available immediately and it can be delivered and installed at your facility with a comprehensive 12 months warranty and extended service coverage

Accessories, Injectors, UPS and Printers are also available upon request. Remember to ask about our PACS systems.